



# LeadCare<sup>®</sup> II

LeadCare is a Magellan Biosciences business

Some healthcare facilities and state policies require users to demonstrate proof of competency on the LeadCare II system. You can print out copies and complete this form to help document each user's training.

**Name:** \_\_\_\_\_ **Facility:** \_\_\_\_\_

## Training Overview

		Date/Initials
<b>LeadCare II Analyzer Kit</b>	<ul style="list-style-type: none"> <li>• <b>Unpack the analyzer kit and review its contents:</b> <ul style="list-style-type: none"> <li>• Read the LeadCare II User's Guide</li> <li>• Read the LeadCare II Quick Reference Guide</li> <li>• Watch the CDC Blood Collection Video and the LeadCare II Instructional Video (approx. 10 minutes each)</li> </ul> </li> </ul>	
<b>LeadCare II Test Kit</b>	<ul style="list-style-type: none"> <li>• Open test kit and review contents</li> <li>• Read the LeadCare II Package Insert</li> </ul>	
<b>Quality Control</b>	<ul style="list-style-type: none"> <li>• Follow the LeadCare II Quick Reference Guide and test a Level 1 and Level 2 Quality Control</li> <li>• Record your results and verify that they are within the acceptable range (provided on the control vial label)</li> </ul>	
<b>User Certification</b>	<ul style="list-style-type: none"> <li>• Review the Skills Checklist</li> <li>• Take the online User Certification Exam</li> </ul>	
<b>Training Completion</b>	<ul style="list-style-type: none"> <li>• Print the training overview &amp; user certificate</li> <li>• Keep for your training documentation files</li> </ul>	

### Need help or additional information?

Contact the LeadCare II Technical Support Team at **(800) 275-0102, ext. 2**, or e-mail **BloodLeadTechSupport@esainc.com**.